

BATH ACCESSORIES



Single Probe Kit.

Mechanical Support for Probes

When setting up a new calibration bath, you need a way to suspend your probes in the bath fluid. We recommend our modular mechanical support systems. Made of only the finest-quality steel and machined parts, these support components combine in hundreds of ways to solve almost any probe suspension problem.

Single Probe Kit

Our single probe kit is a good way to get started. It has one medium clamp, one 10-inch rod, one bosshead, and one V-base for holding one probe. Excellent for Hart's Surface Dry-Well.

Economy Kit

Our economy kit includes two V-bases, one 29-inch stainless steel rod, two 23-inch rods, five bossheads, two micro-clamps, and one medium clamp. This will build a bench-mounted frame for suspending one rod and three clamps over a bath opening. Simply add to your setup as needed. Choose from any of the listed accessories.

Individual Hardware

Our selection of clamps and stands provides a simple way to hold probes and thermometers in baths during calibration. System components can be assembled in a number of ways to suit individual needs. Select clamps, rods, bossheads, and bases to fit your individual needs.

Ordering Information		
Selected Sets		Model Price
Single Probe Kit	Includes: 1 medium clamp, 1 10-inch rod, 1 bosshead, and 1 V-base	#2051
Economy Kit	Includes: 2 V-bases, 1 29-inch stainless steel rod, 2 23-inch rods, 5 bossheads, 2 micro-clamps, and 1 medium clamp	#2050
Individual Hardware		
Micro-Clamp	Holds thermometers and probes with diameters to 0.75 inch.	#2055 (Pkg. of 2)
Medium Clamp	Holds diameters up to 1.75 inches.	#2056 (Pkg. of 1)
Non-Slip Tape	Increases grip on clamps and bossheads.	#2057 (1 roll)
Stainless Steel Rods	Used to assemble supports, frameworks, or scaffolds.	#2058 10" (Pkg. of 1)
		#2059 20" (Pkg. of 1)
		#2060 23" (Pkg. of 1)
		#2061 29" (Pkg. of 1)
Bosshead	Clamps two rods at right angles. Also attaches clamps to rods.	#2062 (Pkg. of 5)
Screw Base	Holds one rod and is screwed to the surface of your bench or bath lid. 2.5-inch-diameter base; screws included.	#2063 (Pkg. of 4)
V-Base	Holds one rod. Weighted for excellent stability. 2.2 pounds (1 kg).	#2064 (Pkg. of 1)
Large V-Base	Same as above but larger. Recommended for holding SPRTs and large probes. 4.4 pounds (2 kg).	#2065 (Pkg. of 1)